

Mindray VS800A Vital Signs Monitor

The VS800A Mindray Vital Signs Monitor measures SpO₂, Temperature, and records data. Suitable for adult, pediatric and neonatal patients. Stable design prevents unit from slipping and tilting back when pressing the buttons.

SmartTemp Thermometer comes with two types of probes: Oral/axillary probe (blue) and rectal probe (red). Probe can auto-warm up to an initial temperature to be sure of an accurate measurement. Neonatal patients are suitable for axillary measurement.

Powerful Data Storage and Management Functions allow you to easily keep and view patient data. Total storage capacity: 100 patient IDs; 1200 sets of data for each parameter. All data can be transferred to a PC through Windows-based software for review and A4 paper printing. Bar code scanner allows you to admit a patient quickly by using a bar case scanner, which is compatible with 1D & 2D bar codes.



- ▶ FDA Approved Masimo SETTM / Nellcor OxiMax™ / Mindray SpO₂
- ▶ SmartTemp™ Thermometer
- ▶ Removable and rechargeable battery, up to 10-hour working time
- ▶ Integrated handle and full mounting solutions
- ▶ Adjustable audible and visual alarms with alarm light indicator
- ▶ High resolution 3.4" LCD and 0.8" large LED displays
- ▶ 3 year warranty
- ▶ 6.7in (D) x 6.7in (W) x 9.4in (H)
- ▶ 7.7 lbs including battery

Item #VSM-800A Vital Signs Monitor with Accessories

Accessories include:

Bar Code Scanner, 2 probes (Adult/Pediatric), 1 NIBP tube (3m) with connector, 1 adult cuff (25cm-35cm), 1 Adult reusable SpO₂ sensor (Mindray), Mindray extension cable (with 6 PIN connector) SmartTemp Acc. Kit, including: 1 Oral/Axillary Probe (blue), 2 boxes Probe covers, Built-in Lithium Ion Battery, & Operation Manual

HeartStart OnSite AED

HeartStart OnSite is an automated external defibrillator created for the everyday person. Available without a prescription, the HeartStart OnSite employs advanced technology designed to make it the easiest to use and most reliable defibrillator.



- ▶ Lightweight: Just 3.3 lbs fully equipped.
- ▶ Easy to Use: Clean design and clear voice instructions, including CPR coaching, instill the confidence that's needed when treating a person in cardiac arrest.
- ▶ Versatile: Available for use on anyone of any age, including children and infants.
- ▶ Advanced Technology: Proven Philips technologies, SMART Analysis and SMART Biphasic, assess heart rhythm and deliver defibrillation.
- ▶ Easy to Maintain: Long life battery and automatic self tests help ensure continued readiness.

Item #AED-1000 HeartStart OnSite AED