

# Anesthesia and Respiratory Carts

## Classic Series Anesthesia Cart 4 Drawer Short or 6 Drawer Tall Workstation Standard Package

All Steel Construction, Key locking system, Vinyl-padded full wrap-around steel bumper, Five-inch (12.7 cm) premium rubber single-wheel caster set, full swivel – one with brake, one directional, Plastic removable top, Aluminum push handles (2), Ball-bearing drawer slides, Double-walled drawer faces with no rough edges, Cart is delivered fully assembled, Certified "Latex Free", Five Year Warranty, Compatible with most accessories, Hard-baked powder coat finish.



### Item #C-6350 4 Drawer Short Workstation

Drawer configuration: 2 x 3-inch (7.6 cm) - 1 x 6-inch (15.2 cm) - 1 x 12-inch (30.5 cm)

Drawer dimensions: 23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled cart dimensions: 35" H x 32" W x 22" D (90.2 cm H x 81.3 cm W x 55.9 cm D)



### Item #C-6450 6 Drawer Tall Workstation

Drawer configuration: 4 x 3-inch (7.6 cm) - 1 x 6-inch (15.2 cm) - 1 x 12-inch (30.5 cm)

Drawer dimensions: 23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled cart dimensions: 41.5" H x 32" W x 22" D (105.4cm H x 55.9cm D x 81.2cm W)

---

## Secure Storage Anesthesia Cart with Breakaway Lock Standard Package



All Steel Construction, Tall cabinet with stainless steel breakaway locking system, Breakaway lock seals, pkg-50, Vinyl-padded full wrap-around steel bumper, Five-inch (12.7 cm) sealed ball-bearing premium rubber single-wheel casters, full swivel – one with brake, one directional, Aluminum push handle/side rail, Removable plastic top, Ball-bearing full extension drawer slides, Double-walled drawer faces with no rough edges, Cart is delivered with wheels installed; some accessories require user installation. Certified "Latex Free", Five Year Warranty, Compatible with most accessories, Hard-baked powder coat finish.

### Item #C-6600 Secure Storage Anesthesia Cart with Breakaway Lock

Drawer configuration: 6 x 3-inch (7.6 cm) - 2 x 6-inch (15.2 cm)

Drawer dimensions: 23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled cart dimensions: 41.5" H x 32" W x 22" D (105 cm H x 81.3 cm W x 55.9 cm D)

Twelve available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD),  
Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and  
Hammer Tone Green (-HTGN)

# Deluxe Anesthesia, Mini Line Anesthesia, & MR Safe Cart

## Deluxe Anesthesia Cart



All Steel Construction, Raised back rail system, 6-Compartment tilt bin organizer, Tape and label dispenser, Large utility basket, Catheter box, 6-Socket hospital grade outlet with 15-foot cord and cord wrap, Plastic waste container, Large sharps container and bracket, 3-inch Exchange tray (2), 6-inch Exchange tray (2), 12-inch Adjustable drawer divider set, Tall cabinet with key locking system, Vinyl-padded full wrap-around steel bumper, Five-inch (12.7 cm) premium rubber single-wheel caster set, full swivel – one with brake, one directional, Plastic removable top, Aluminum push handles (4), Ball-bearing drawer slides, Double-walled drawer faces with no rough edges, Cart is delivered wheels assembled, Certified “Latex Free”, Five Year Warranty, Compatible with most accessories, Hard-baked powder coat finish

Item #C-6456 Deluxe Anesthesia Cart Classic Line

Drawer configuration: 4 x 3-inch (7.6 cm) - 1 x 6-inch (15.2 cm) - 1x12-inch (30.5 cm)

Drawer dimensions: 23" W x 17" D (58.4 cm x 43.2 cm)

Assembled cart dimensions: 41.5" H x 32" W x 22" D (105 cm H x 81.3 cm W x 55.9 cm D)

## Mini Line Anesthesia Cart Standard Package



Painted Steel Cabinet and Drawers, Space Saver Design fits into smaller areas, Allows supply storage in areas where space is limited, Cabinet with 30 vertical inches (76.2 cm) drawer space, Keylocking system with two keys, 3-inch (7.6 cm) single-wheel stem caster set, full swivel – two with brake, Pre-threaded accessory mounting holes, Ball-Bearing Full extension Drawer Slides, Removable Plastic Top, Certified “Latex Free”, Five Year Warranty, Cart requires caster installation and accessory assembly after delivery, Hard-baked powder coat finish.

Item #C-3156 Mini Line Anesthesia Cart

**Drawer configuration:** 3 x 3-inch (7.6 cm) - 2 x 6 -inch (15.2 cm) - 1 x 9-inch (22.9 cm)

**Assembled Cart Dimensions:** 38.25" H x 18" W x 18" D (97.2 cm H x 45.7cm W x 45.7 cm D)

## MR Safe Cart



All-Aluminum Cabinet and Drawers, Tall Cabinet (33.25 inches / 84.5 cm) with Breakaway lock, Vinyl-padded full wrap-around steel bumper, Five-inch (12.7 cm) MR-Safe single-wheel caster set, full swivel – two with brake, Stainless Steel Pull-out Shelf, Pre-threaded accessory mounting holes, Low Ferrous Ball-Bearing Drawer Slides, Replaceable Plastic Top with molded-in Push Handles, Cart is delivered fully assembled.

Certified “Latex Free”, Certified MR-Safe by independent testing authority (copies of report available), Five Year Warranty, MR-Safe accessories available, Hard-baked powder coat finish. Shelf Painted Gloss White with MRI Sticker.

Item #C-MR6B Drawer MR-Safe Anesthesia Cart with Breakaway Lock


Drawer configuration: 1 x 4-inch (10.2 cm) - 2 x 3.25-inch (8.3 cm) - 2 x 6.5-inch (16.5 cm) - 1 x 9.75-inch (24.8 cm)

Assembled Cart Dimensions: 45.7 5" H x 30.5" W x 23.75" D (116.2 cm H x 77.5 cm W x 60.3 cm D)

# Accessories Options and Replacement Parts



C680302



C680431



C680404




C680419



C680409



C68157



C680415



C40108

## Accessories Options & Replacement Parts for Anesthesia Carts:

- C680404 Electrical Outlet, Hospital Grade, 6 Socket, 15 Foot Cord, Includes Cord Wrap
- C68161 Aluminum Push Handle, 16"
- C68164 Back Rail For Monitor Cart Shelf Support
- C68157 Back Rail System With 3-24" Rails, Raised, Classic Line
- C68163 Back Rail System With 2-24" Rails, Raised, Classic Line
- C40088 Breakaway Seals Unnumbered pk/50
- C680311 Short Overhead Shelf Support Set
- C680301 Side Mounted Drop Shelf With Lip
- C680425OHMS IV Pole - 2 prong For use with Short Overhead Monitor Shelf OHMS
- C680419 Rail Clip
- C680415 Tape/Label Dispenser
- C683803 Tilt Bin, 4 Compartments, With Rail Clips
- C683804 Tilt Bin, 5 Compartments, With Rail Clips
- C683805 Tilt Bin, 6 Compartments, With Rail Clips
- C683802 Sharps Container And Mounting Bracket, 3 gallon
- C608401 IV Pole - 2 prong
- C680400 IV Pole - 4 prong
- C680409 Steel waster container with cover
- 680431 Utility basket with brackets

## Accessories Options & Replacement Parts for Mini Line Anesthesia Cart:

- CMINI-RBR Back Rail System with 2-18" Rails, Raised, Mini Line
- C680404 Electrical Outlet, Hospital Grade, 6 Socket, 15 Foot Cord, Includes Cord Wrap
- C68161 Aluminum Push Handle, 16"
- C680408 Waste Container Without Cover, Plastic, 8 Qt.
- C683803 Tilt Bin, 4 Compartments, With Rail Clips
- C683804 Tilt Bin, 5 Compartments, With Rail Clips
- C683805 Tilt Bin, 6 Compartments, With Rail Clips
- C68620-P1 Mini-Line Drawer Divider Tray Pack 1 Consisting Of 3 Long, 3 Short.  
1 (1/2)L, 1 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- C68620-P2 Mini-Line Drawer Divider Tray Pack 2 Consisting Of 5 Long, 5 Short.  
2 (1/2)L, 2 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- C683802 Sharps Container And Mounting Bracket, 3 gallon
- CMINI-DRW9N Mini Line 9" Drawer With Locking Narcotics Compartment

## Accessories Options & Replacement Parts for MR-Safe Anesthesia Cart:

- CMR-RBR2 2 Tier Raised Back Rail, MR-Safe Cart
- CMR-SR Side Rail, MR-Safe Cart
- CMR-CARDBRD Cardiac Board And Mounting Brackets, MR-Safe Cart
- CMR-RLHK Rail Hook, MR-Safe
- CMR-IV2 IV Pole, MR-Safe
- CMR-O2 O2 Bracket, MR-Safe
- CMR-UHOOKS4 Utility Hooks For MR-Safe Carts (Set Of 4)
- CMR-4TILT 4-Compartment Tilt Bin, MR-Safe, with mounting clip
- CMR-5TILT 5-Compartment Tilt Bin, MR-Safe, with mounting clip
- CMR-6TILT 6-Compartment Tilt Bin, MR-Safe, with mounting clip
- CMR-3GWASTE Waste Container, 3 Gallon, with mounting bracket, MR-Safe



C680310



C683802



C68163

Please call Customer Service 1-888-883-7331 for a complete listing of Accessories Options & Replacement Part Accessories.