

Cricothyrotomy Kits

Quick Trach II with Cuff

VBM

Cricothyrotomy is an emergency procedure to avoid death from suffocation in case of upper airway obstruction if intubation or tracheotomy are impossible (e.g. glottis edema, foreign body, carcinoma of the larynx).



Product Specifications:

- ▶ Size 4.0 mm
- ▶ Preassembled and ready to use
- ▶ Sterile, Single Use, Latex Free

Product Benefits:

Safety Clip

- ▶ The plastic cannula is pushed forward until the safety clip clicks into position
- ▶ The safety clip is firmly fixed at the connector and therefore avoids that the metal needle is pushed out of the cannula again by mistake
- ▶ The needle tip is covered by the plastic cannula and therefore does not cause injury to the posterior tracheal wall

Necktape

- ▶ For soft fixation
- ▶ Soft Foam material

Stopper

- ▶ Prevents the needle from being inserted too deep and therefore reduces the risk of posterior tracheal wall perforation

Cuff

- ▶ Allows sufficient ventilation with aspiration protection
- ▶ Cuff is made of ultra-thin material and is very durable

Components:

- ▶ 1 Quick Trach, sterile with stopper
- ▶ 1 Syringe
- ▶ 1 Cuff
- ▶ 1 Cushion neckband

1-30104 Quick Trach II with Cuff Emergency Cricothyrotomy Kit 4.0mm

Quick Trach

Quick Trach allows quick and safe access for ventilation in the presence of acute respiratory distress with upper airway obstruction.



Product Specifications:

- ▶ Size 1.5mm, 2.0mm or 4.0 mm
- ▶ Plastic cannula with fixation flange and 15 mm connector
- ▶ Removable stainless steel cricothyrotomy needle
- ▶ Removable stopper
- ▶ Preassembled and ready to use
- ▶ Sterile
- ▶ Single Use
- ▶ Latex-free

RUSCH

Product Benefits:

- ▶ Removable stopper reduces risk of damage to the rear wall of the trachea
- ▶ Conical needle tip guarantees the smallest necessary stoma and reduces bleeding
- ▶ The syringe allows identification of the trachea by aspirating air
- ▶ Connection tube for immediate ventilation

Components:

- ▶ 1 Quick Trach, sterile with stopper
- ▶ 1 Syringe
- ▶ 1 Connecting tube with 15mm adapter
- ▶ 1 Cushion neckband

Baby 1-20905 Quick Trach Emergency Cricothyrotomy Kit 1.5mm
Child 1-20902 Quick Trach Emergency Cricothyrotomy Kit 2.0mm
Adult 1-2090 Quick Trach Emergency Cricothyrotomy Kit 4.0mm