

Cuff Pressure Accessories

Cuff Pressure Gauge

VBM

Pocket



- ▶ Inflates and controls cuff pressure for tracheal tubes
- ▶ Use with high volume, low pressure cuffs
- ▶ Scale reads 0-120cm H2O
- ▶ Protection ring included
- ▶ Supplied with 1m (40")' connection tubing

Item #7-5400

Pediatric with Hook



- ▶ Inflates and controls cuff pressure for tracheal tubes
- ▶ Release valve allows sensitive adjustment of the desired cuff pressure
- ▶ 50mm scale
- ▶ Pressure range from 0-60cm H2O
- ▶ Protection ring included
- ▶ Supplied with 1m (40")' connection tubing
- ▶ Includes hook

Item #7-5401

Variable with Adjustable Pointer



Silicone tracheal tubes have a high permeability of anesthetic gases. The minimal necessary cuff pressure can be marked with the adjustable pointer. Pressure increase can be released to the marked pressure with the pressure release valve.

- ▶ Inflates and controls cuff pressure for tracheal tubes
- ▶ Use with high volume, low pressure cuffs
- ▶ Scale reads 0-120cm H2O
- ▶ Protection ring included
- ▶ Supplied with 1m (40")' connection tubing

Item #7-5460

Cuff Pressure Monitors



- ▶ Gauge & Adapter
- ▶ Reusable

Item #7-NS-CPR
Item #7-NS-CPB

+/- 60cm H2O rear mount 1/4" NPT
+/- 60cm H2O bottom mount 1/4" NPT