

# Cylinder Carts and Racks



- ▶ Manufactured from anodized gray structural tubing
- ▶ 6 1/2" wheels
- ▶ Lift handle
- ▶ Adjustable height
- ▶ Grip handle
- ▶ Hook for hose or mask
- ▶ Bolt to secure cylinder
- ▶ Transports D & E cylinders

Item #7001



- ▶ Manufactured from durable structural tubing
- ▶ Holds 2 cylinders
- ▶ 6 1/2" wheels
- ▶ Durable powder coat finish
- ▶ Grip handle
- ▶ Hook for tubing or mask
- ▶ Bolt to secure cylinders
- ▶ Transports D & E cylinders

Item #7002



- ▶ Fabricated from structural tubing
- ▶ Heavy Duty wheels & Axle
- ▶ Durable powder coat finish
- ▶ Safety chain included

Item #7005  
(Transports M, M60 & H)

Item #7005M  
(Transports M & M60)

Also Available with Retractable and Non-Retractable Wheels



- ▶ Transport or stationary rack for M7, M9, C, D & E cylinders
- ▶ Manufactured from structural tubing
- ▶ Durable finish
- ▶ Heavy Duty wheels & axle


Item #7012  
Holds 12 cylinders

Item #7012-6  
Holds 6 cylinders

Item #7012-4  
Holds 4 cylinders

Item #7012-W  
Holds 12 cylinder w/ 4 wheels

## Cylinder Racks



Capacity	Cylinder	Part #
4	C, D, or E	7004
6	C, D, or E	7006
6	M7, M9, C, D, & E	7006-M
12	M6/B	7029
12	C, D, or E	7010
12	M7, M9, C, D, & E	7010-M
18	M7, M9 C, D, or E	7003-M
24	M7, M9, C, D, or E	7008
40	D, or E	7009

Wall Mount  
Holds 1 D or E Cylinder  
Item #7000



Are you looking for a rack or cart not shown?  
We offer many more options to choose from.  
Inquire today, 1-888-883-7331.

# MRI Compatible Products

## Flowmeter Features:


- Non-Magnetic, Slim Line model
- Anodized aluminum body with brass adjustment knob
- Polycarbonate inner and outer tubes
- Calibrated at 50 psi inlet
- DISS 1240 outlet
- 6 year limited warranty

## Regulator Features:

- Non-Magnetic Preset and Flowmeter Regulator models available
- All brass body and high pressure chamber
- Internal reseating relief valve protects against over pressurization
- 2" diameter gauge
- Sintered filter for additional safety and to extend regulator life
- Durable neoprene diaphragm
- Maximum inlet pressure 3000 psi
- DISS 1240 outlet
- 6 year limited warranty



Slim Line Flowmeter  
Item #7501-MRI



Preset Pressure Gauge Regulator  
50 psi Preset, CGA-870 Yoke  
Item #7770-MRI



Single - Stage Flowmeter Regulator  
½ - 15 LPM, 50 psi Present, CGA-870 Yoke  
Item #7715-MRI

## MRI Compatible Carts & Rack (Non-Magnetic)



Item #7004-MRI  
Holds 1 (H or T) Cyl.



Item #7005-MRI  
Holds 1 (M or M60) Cyl.



Item #7001-MRI  
Holds 1 (D or E) Cyl.



Item #7002-MRI  
Holds 2 (D or E) Cyl's



Item #7000-MRI  
Wall Mount  
Holds (D or E) Cyl's