


# Dental Guard, Fiber Optic German Profile Laryngoscopes, Disposable Laryngoscopes & Fiber Optic Stylettes

## Dental Guard by CCR Medical

- ▶ Protects against damage to teeth during intubation, bronchoscopy, laryngoscopy, tonsillectomy, and tooth extraction
- ▶ Accommodates adults & children
- ▶ Blue, Non-toxic vinyl plastic
- ▶ Latex Free

LATEX FREE




Item #1990 Dental Guard  
Item #1990-10 Dental Guard

REQUEST A SAMPLE



## Malleable Fiber Optic Stylette

- ▶ Aids in the proper placement of E/T tubes
- ▶ Incorporates the features of a malleable stylette
- ▶ Illuminates the soft tissue of the neck
- ▶ Uses 2 "AA" batteries



3675-14 Pediatric 14" Fiberoptic Stylette  
3675-18 Adult 18" Fiberoptic Stylette  
1007 2.5V .35mm Replacement Bulb

## German Profile Fiber Optic Laryngoscopes



Large oval integrated bundle for optimum light configuration and output  
303/304 surgical stainless steel  
Interchangeable with all Green Fiber Optic Systems  
Up to 50% more light

### MacIntosh Profile

2195-0	Newborn	Size 0
2195-1	Infant	Size 1
2195-2	Child	Size 2
2195-3	Adult	Size 3
2195-4	Large Adult	Size 4

### Miller Profile

2295-00	Premature	Size 00
2295-0	Newborn	Size 0
2295-1	Infant	Size 1
2295-2	Child	Size 2
2295-3	Adult	Size 3
2295-4	Large Adult	Size 4

**New Product**

## Disposable Metal Fiberoptic Laryngoscope Blades


Metal-Disposable for single use. Non-Sterile. The disposable laryngoscope blade that fits on all green color coded fiberoptic handles. Bright white light out-performs other blades. Tested for tensile strength.

### MacIntosh Profile



2-5512	Size 2	100mm
2-5513	Size 3	130mm
2-5514	Size 4	155mm

### Miller Profile



2-5500	Size 0	75mm
2-5501	Size 1	102mm
2-5502	Size 2	155mm