

ECG Electrodes



Cleartrace® - Adult Electrode



Clear NeoDerm® tape allows visual monitoring of the patient's skin condition at all times. Perforated tape permits skin ventilation and moisture evaporation for increased patient comfort. Conductive adhesive gel. Conforms to patients body without pulling on the skin. Will not absorb fluids. Latex free.

11-1703 3/pouch
11-1705 5/pouch

Cleartrace® 2 - Adult Electrode MRI Compatible



Square/Oval shape adds "staying power" to the design - up to 7 days of use. Clear NeoDerm® tape allows visual monitoring of the patient's skin condition at all times. NeoDerm® tape is waterproof, yet allows the skin to breath. Available with a carbon snap for radiotranslucent monitoring. Latex free.

11-2703 3/pouch
11-2705 5/pouch

Pediatric/Neonatal Tape Snap Electrodes



Clear NeoDerm tape allows visual monitoring of the baby's skin condition at all times. Low chloride gel and light adhesive coat weight reduces the possibility of skin irritation. Perforated tape permits skin ventilation and moisture evaporation to increase comfort. Tab provides for easy application and removal. Wet get formulation is gentle to the skin. Latex free.

11-1601 1/pouch
11-1603 3/pouch

Softrace® Pediatric and Neonatal Electrodes Wire Attached ECG Electrodes



Softrace® soft fabric backing and mild conductive adhesive provide secure, comfortable electrode contact. Small size and extremely low profile increases your placement options. Humidity has little effect on Softrace® electrodes' ability to stick or to be removed without damaging fragile skin. Latex Free.

11-2321 Softrace® Neonatal 3/pouch
11-2331 Softrace® Pediatric 3/pouch
11-2310 Softrace® Limb Band 3/pouch



PadPro Adult and Pediatric Electrodes

Defibrillation - Pace - Cardioversion - Monitoring



Adult Radiotransparent

Non-Sterile, Latex Free, 13.02 X 8.57 cm

11-2001Z	Zoll
11-2001ZPC	Zoll Pre-Connect
11-2001ZC	Zoll Pad and Wires
11-2001M	Medtronic
11-2001MPC	Medtronic Pre-Connect
11-2001MC	Medtronic Pad and Wires
11-2001P	Philips
11-2001PPC	Philips Pre-Connect
11-2001PC	Philips Pad and Wires

Adult Radiotranslucent

Non-Sterile, Latex Free, 13.02 X 8.57 cm

11-2516Z	Zoll
11-2516ZPC	Zoll Pre-Connect
11-2516M	Medtronic
11-2561MPC	Medtronic Pre-Connect
11-2516P	Philips
11-2516PPC	Philips Pre-Connect

Pediatric Radiotranslucent

Non-Sterile, Latex Free, 7.6 X 8.2 cm

11-2603Z	Zoll
11-2603M	Medtronic
11-2603P	Philips



Mini-Pediatric Radiotranslucent

Non-Sterile, Latex Free, 5.85 X 5.85 cm

11-2602Z	Zoll
11-2602M	Medtronic
11-2602P	Philips

