


Portable Emergency Oxygen Manifold Systems

OCTI-FLO2 Portable Emergency Oxygen Manifold




- ▶ Use for mass casualty, EMS, and triage situations
- ▶ Oxygen delivery for up to 8 patients
- ▶ Flowmeters are individually controlled
- ▶ 8 easy-to-read click-style flowmeters from 0-25 LPM
- ▶ Pressure gauge on manifold
- ▶ Standard DISS 1240 inlet and a 100 psi relief valve safety feature
- ▶ Lightweight and compact durable aluminum construction
- ▶ Includes cannulas, breathing masks, and compact carry case
- ▶ Affordable high quality craftsmanship for years of dependable service



Item #10-EM08 Emergency Oxygen Manifold
Complete Kit

Includes: Cannulas, Breathing Masks, and Compact Carry Case

Optional Components to Complete the Emergency Oxygen Delivery System



(0-100 psi)

Item #7021-100
0-100 psi Regulator CGA 540 H Tank




(Oxygen medical hose)

Item #HA-C4OXY240FNNFNN
20 foot oxygen medical hose with DISS Fittings

Oxygen Supply Manifold

- ▶ Use for mass casualty, EMS, and triage situations
- ▶ Oxygen delivery for up to 6 patients
- ▶ Flowmeters are individually controlled
- ▶ 8 easy-to-read click-style flowmeters from 0-15 LPM
- ▶ Pressure gauge on manifold
- ▶ Standard DISS 1240 inlet and a 50 psi supply pressure requirement
- ▶ Complete kit weights 6.5 lbs
- ▶ Includes 20 ft. Oxygen Hose, 6 oxygen therapy Masks, and carry case
- ▶ Affordable high quality craftsmanship for years of dependable service

PRECISION MEDICAL



Item #10-EM06 Oxygen Supply Manifold Complete Kit

Includes: 20 ft. Oxygen Hose, 6 Oxygen Therapy Masks, and Carry Case

Item #10-EM68 Oxygen Supply Manifold Only

Optional Component to Complete the Emergency Oxygen Delivery System



Item #10-540
0-50 psi Preset Regulator
CGA 540 H Tank