

# MRI Compatible Products

## Flowmeter Features:

- Non-Magnetic, Slim Line model
- Anodized aluminum body with brass adjustment knob
- Polycarbonate inner and outer tubes
- Calibrated at 50 psi inlet
- DISS 1240 outlet
- 6 year limited warranty

## Regulator Features:

- Non-Magnetic Preset and Flowmeter Regulator models available
- All brass body and high pressure chamber
- Internal reseating relief valve protects against over pressurization
- 2" diameter gauge
- Sintered filter for additional safety and to extend regulator life
- Durable neoprene diaphragm
- Maximum inlet pressure 3000 psi
- DISS 1240 outlet
- 6 year limited warranty



Slim Line Flowmeter  
Item #7501-MRI



Preset Pressure Gauge Regulator  
50 psi Preset, CGA-870 Yoke  
Item #7770-MRI



Single - Stage Flowmeter Regulator  
½ - 15 LPM, 50 psi Present, CGA-870 Yoke  
Item #7715-MRI

## MRI Compatible Carts & Rack (Non-Magnetic)



Item #7004-MRI  
Holds 1 (H or T) Cyl.



Item #7005-MRI  
Holds 1 (M or M60) Cyl.



Item #7001-MRI  
Holds 1 (D or E) Cyl.



Item #7002-MRI  
Holds 2 (D or E) Cyl's



Item #7000-MRI  
Wall Mount  
Holds (D or E) Cyl's