

# Infusion Cuffs and Accessories

## Shock Resistant Gauge - Durable Design



### Reusable Pressure Infusion Cuff 500ml & 1000ml


- ▶ Transparent cuff allows clear view of the fluid bag and the fluid level
- ▶ Wrap-around design assures consistent pressure on the complete cuff
- ▶ Shock resistant precision manometer is over-pressure proved and assures exact pressure control
- ▶ Stop cock allows quick and complete deflation of the cuff
- ▶ Large hand inflator for quick inflation of the cuff
- ▶ Easily clean with mild soap and water

### Reusable Infusion Cuff with Bulb

- 7-57850 Pressure Infusion Cuff 500ml with manometer and hand inflator
- 7-57800 Pressure Infusion Cuff 1000ml with manometer and hand inflator

### Disposable Pressure Infusion Cuff 500ml & 1000ml

- ▶ Durable material (pressures up to 300mmHg)
- ▶ Adjustable hook for bag attachment
- ▶ Transparent netting allows clear view of fluid level
- ▶ Single patient use to prevent cross-contamination
- ▶ Latex Free



### Disposable Infusion Cuff with Bulb

#### 500 ml Cuffs

- 7-56205 Pressure Infusion Cuff with (B) hand inflator, release valve, & color coded gauge
- 7-56905 Pressure Infusion Cuff with (C) \*manometer gauge, bulb, & valve
- 7-56305 Pressure Infusion Cuff with (D) hand inflator and release valve

#### 1000 ml Cuffs

- 7-56201 Pressure Infusion Cuff with (B) hand inflator, release valve, & color coded gauge
- 7-56901 Pressure Infusion Cuff with (C) \*manometer gauge, bulb, & valve
- 7-56301 Pressure Infusion Cuff with (D) hand inflator and release valve

\*reusable



### Accessories for Infusion Cuff

- 7-5103 Universal Adjustable 30" Telescopic Infusion Pole with Clamp

### Spare Accessories and Parts for Reusable Infusion Cuff

- 7-3470 (E) Manometer 300mmHg.
- 7-5136 (D) Hand Inflator with release valve and male luer lock

### Spare Accessories and Parts for Disposable Infusion Cuff

- 7-5610 (A) Pressure Infusion Cuff 500ml with female luer lock
- 7-5611 (A) Pressure Infusion Cuff 1000ml with female luer lock
- 7-5126 (B) Hand Inflator with release valve, color coded gauge and male luer lock
- 7-5196 (C) Manometer 300mmHg, with hand inflator, release valve and male luer lock
- 7-5136 (D) Hand Inflator with release valve and male luer lock