


CCR Medical Fiberoptic Green System Laryngoscopes & Handles

- ▶ Supplied with a halogen bright bulb
- ▶ Chrome plated brass or stainless steel
- ▶ Steam autoclaveable or cold sterilization
- ▶ Available in 4 popular sizes
- ▶ Completely interchangeable with all green system and fiberoptic blades and handles

- 2750-1 Penlight Fiberoptic C/P "AA"
- 2750-2 Stubby (Bantam) Fiberoptic C/P "AA"
- 2750-3 Medium Fiberoptic C/P "C"
- 2750-4 Large Fiberoptic C/P "D"
- 2750-5 Howland Lock

1990 Halogen Fiberoptic Bulb 2.5 volt



Howland Lock

The Howland Lock is designed to assist in the angle placement of the blade for bull-neck and obese patients. Recommend use with straight blades.



2750-5



German Profile Fiber Optic Laryngoscopes

Large oval integrated bundle for optimum light configuration and output
303/304 surgical stainless steel
Interchangeable with all Green Fiber Optic Systems
Up to 50% more light

MacIntosh Profile

2195-0	Newborn	Size 0
2195-1	Infant	Size 1
2195-2	Child	Size 2
2195-3	Adult	Size 3
2195-4	Large Adult	Size 4

Miller Profile

2295-00	Premature	Size 00
2295-0	Newborn	Size 0
2295-1	Infant	Size 1
2295-2	Child	Size 2
2295-3	Adult	Size 3
2295-4	Large Adult	Size 4

- ▶ Manufactured from 303/304 surgical stainless steel
- ▶ Interchangeable with all green F/O systems
- ▶ Steam autoclaveable or cold sterilization
- ▶ 3 millimeter fiber bundle for optimum light out put
- ▶ Halogen bright light

LATEX FREE



MacIntosh American Profile



2050-0	Premature	Size 0	80mm
2050-1	Infant	Size 1	89mm
2050-2	Child	Size 2	110mm
2050-3	Adult	Size 3	126mm
2050-35	Adult	Size 3.5	144mm
2050-4	Large Adult	Size 4	162mm

MacIntosh English Profile Fiberoptic



2150-1	Infant	Size 1	92mm
2150-2	Child	Size 2	100mm
2150-3	Adult	Size 3	130mm
2150-4	Large Adult	Size 4	155mm

Miller Profile Fiberoptic



2250-0	Newborn	Size 0	78mm
2250-1	Infant	Size 1	105mm
2250-15	Small Child	Size 1.5	128mm
2250-2	Child	Size 2	158mm
2250-3	Medium Adult	Size 3	198mm
2250-4	Adult	Size 4	208mm

Phillips™ American Profile Fiberoptic



2450-1	Child	Size 1	113mm
2450-2	Adult	Size 2	166mm

Robertshaw English Profile Fiberoptic



2550-0	Infant	Size 0	94mm
2550-1	Child	Size 1	112mm

Wis-Hipple/Wisconsin American Profile Fiberoptic




2650-0	Newborn	Size 0	76mm
2650-1	Infant	Size 1	102mm
2650-15	Small Child	Size 1.5	115mm
2650-2	Child	Size 2	135mm
2650-3	Adult	Size 3	162mm
2650-4	Large Adult	Size 4	199mm

CCR Medical Conventional Laryngoscopes & Handles

COMPARE OUR QUALITY AND PRICE!

- ▶ Compatible with all other brands
- ▶ Chrome plated brass or stainless steel
- ▶ Steam autoclaveable or cold sterilization
- ▶ Available in 4 popular sizes
- ▶ Uses "AA", "C", or "D" cell batteries

- 2700-1 Penlight C/P Conventional "AA"
- 2700-2 Stubby (Bantam) C/P Conventional
- 2700-3 Medium C/P Conventional "C"
- 2700-4 Large C/P Conventional "D"
- 2700-5 Howland Lock



Howland Lock
The Howland Lock is designed to assist in the angle placement of the blade for bull-neck and obese patients. Recommend use with straight blades.



2700-5

MacIntosh American Profile



2000-0	Premature	Size 0	80mm
2000-1	Infant	Size 1	89mm
2000-2	Child	Size 2	108mm
2000-3	Adult	Size 3	128mm
2000-35	Adult	Size 3.5	144mm
2000-4	Large Adult	Size 4	159mm

MacIntosh English Profile



2100-0	Premature	Size 0	83mm
2100-1	Infant	Size 1	92mm
2100-2	Child	Size 2	100mm
2100-3	Adult	Size 3	130mm
2100-4	Large Adult	Size 4	155mm

Miller American Profile



2200-00	Premature	Size 00	65mm
2200-0	Newborn	Size 0	77mm
2200-1	Infant	Size 1	102mm
2200-15	Small Child	Size 1.5	126mm
2200-2	Child	Size 2	155mm
2200-3	Adult	Size 3	195mm
2200-4	Large Adult	Size 4	205mm

Miller English Profile



2300-0	Newborn	Size 0	77mm
2300-1	Infant	Size 1	102mm
2300-2	Child	Size 2	153mm
2300-3	Adult	Size 3	195mm
2300-4	Large Adult	Size 4	205mm

- ▶ Compatible with all other brands
- ▶ Five year warranty (excluding bulb)
- ▶ Manufactured from 303/304 surgical stainless steel
- ▶ Steam autoclaveable or cold sterilization
- ▶ Superior quality

LATEX FREE

Phillips™ American Profile



2400-1	Child	Size 1	113mm
2400-2	Adult	Size 2	166mm

Robertshaw English Profile



2500-0	Infant	Size 0	94mm
2500-1	Child	Size 1	112mm

Wis-Hipple American Profile




2600-0	Premature	Size 0	76mm
2600-1	Infant	Size 1	102mm
2600-15	Small Child	Size 1.5	115mm

Wisconsin American Profile



2600-2	Child	Size 2	135mm
2600-3	Adult	Size 3	162mm
2600-4	Large Adult	Size 4	199mm

- 1001 Small Frosted
- 1002 Large Frosted
- 1003 Small Clear
- 1004 Large Clear



Laryngoscope Repair Service

Conventional Blades: Conventional Handles:

Mail to: CCR Medical - Attn: Repair - 967 43 Avenue NE - St. Petersburg, FL 33703