


Lighted Stylette and Disposable Metal Laryngoscope Blades

The Fiber Optic Malleable Lighted Stylette provides illumination during intubation and also incorporates the features of a malleable stylette. Assists in the visualization required for the proper placement of tracheal tube during intubation and effectively illuminates the soft tissue of the neck.

Lighted Stylette

- 3675-14 Pediatric Stylette with Handle
- 3675-18 Adult Stylette with Handle
- 3675-S Stylette Stop, delrin
- 3675-14R Pediatric 25cm Replacement Stylette
- 3675-18R Adult 36cm Replacement Stylette
- 3675-H Replacement Handle 2 "AA" Batteries
- 3675-B Replacement Bulb Extra Bright



Disposable Metal Conventional Laryngoscope Blades

Metal-Disposable for single use. Non-Sterile. Individually packaged.

The disposable laryngoscope blade is compatible with all conventional handles.

MacIntosh Profile



- 2-5532 Size 2
- 2-5533 Size 3
- 2-5534 Size 4

Miller Profile




- 2-5520 Size 0
- 2-5521 Size 1
- 2-5522 Size 2
- 2-5523 Size 3

Disposable All Metal Fiber Optic Laryngoscope Blades



MacIntosh Profile



- 12-4712 Size 1
- 12-4712 Size 2
- 12-4713 Size 3
- 12-4714 Size 4

Miller Profile



- 12-4500 Size 00
- 12-4520 Size 0
- 12-4521 Size 1
- 12-4522 Size 2
- 12-4523 Size 3
- 12-4524 Size 4

The Venticare range of fiber optic laryngoscopes are all metal and lightweight.

The halogen bulb produces a cold yet bright illumination, providing an intensity of light required during examination and intubation.

Individually package pre-sterilized blade.

Green RuschLite Disposable Metal Fiberoptic Laryngoscope Blades

Metal-Disposable for single use. Non-Sterile. The disposable laryngoscope blade that fits on all green color coded fiberoptic handles. Bright white light out-performs other blades. Tested for tensile strength.




MacIntosh Profile



- 2-5512 Size 2 100mm
- 2-5513 Size 3 130mm
- 2-5514 Size 4 155mm

Miller Profile



- 2-5500 Size 0 75mm
- 2-5501 Size 1 102mm
- 2-5502 Size 2 155mm