

Mask Harness & Reusable Cushion Face Masks


Mask Harness

LATEX FREE

- ▶ Manufactured from 1/16" medical grade latex free nitrile
- ▶ Material is sturdy but designed to fit comfortably
- ▶ Provides support of face mask with complete adjustability
- ▶ White, clean user friendly color
- ▶ Available in Adult, Child, and Infant Sizes
- ▶ Reusable
- ▶ Substantial Cost Savings! Compare.
- ▶ **Made in the USA**

Item #5995 Adult Size
 Item #5996 Child Size
 Item #5997 Infant Size

One of Each Size Item #5999




Bi-Mask Reusable Air Cushion Face Mask

LATEX FREE

This flexible mask and unique design is soft and comfortable for the patient. With a combination of two materials, the plastic shell is molded into the silicone and gives the mask stability allowing the user to manipulate the mask for a perfect seal.

Latex Free - Autoclavable up to 134°C



Without Hook Ring

6-6523	Adult Small	Yellow
6-6524	Adult Medium	Red
6-6525	Adult Large	Purple
6-6526	Adult Ex. Large	Blue

With Hook Ring

6-6523HR	Adult Small	Yellow
6-6524HR	Adult Medium	Red
6-6525HR	Adult Large	Purple
6-6526HR	Adult Ex. Large	Blue

Silicone Air Cushion Face Mask

LATEX FREE

- ▶ Polysulfone hook ring
- ▶ Reusable, Made of Silicone
- ▶ Inflatable
- ▶ Anatomically Shaped Dome
- ▶ 22mm ID connector
- ▶ Autoclaveable up to 134°C

6-3522	Child with Ring	Green	Size 2
6-3523	Small Adult with Ring	Yellow	Size 3
6-3524	Medium Adult with Ring	Red	Size 4
6-3525	Large Adult with Ring	Purple	Size 5
6-3526	Ex Large Adult with Ring	Blue	Size 6

