

# Timeter® PCS 414 Portable Air Compressor



Timeter® PCS 414 Portable Air Compressor is so much quieter than other air compressors that you may wonder if it's running. You'll like its looks, too. The durable, mar-resistant steel case is easy to clean and maintain. The unit weights only 26.5 pounds. Engineered for maximum efficiency and long service life, with an output rating of 14-lpm @ 50 psig. Permanent lubricated heavy-duty bearings are designed to operate 24 hours per days, 365 days per year. Cool-running and lightweight, the PCS 414 is ideal for use with high-pressure and low-pressure aerosol administration devices such as pediatric tents, croupettes, mist tents, and large and small nebulizers, including those being prescribed for the aerosolization of Pentamidine. It delivers safer, continuous high performance, day in an day out.



10-T4614	PCS 414 Air Compressor
10-T6159	Portable Transfer Cart
10-TZ85	Right Angle Adaptor
10-TZ111	Short "J" Bracket
10-TZ120	8" J Bracket
10-TZ70	Cooling Coil
10-TZ72	Condensation Trap

115 volt, 60 Hz	Output: 14 lpm @ 50 psig / 19 lpm @ 20 psig
Dew Point Depression: Ambient Air	Power: 4.6 Amps (nominal)
Power Cord: #18-3 X 6'7"	Plug: Molded Hospital Grade
Electrical Protection: 8.0 Amp Fuses (2)	Ground Impedance: 0.1 Ohms Maximum
Dimensions: 8 1/4" W x 21" D x 13 1/2" H	Weight: 25.6 lbs. (11.61 kg)
Outlet: Oxygen Male DISS	Nose Level: 50 dBa Pressure Gauge: 0-100 psi
Pressure Regulator: 0-50 psi	Approvals: UL and CSA Warranty: 1 year

## Portable Mist Tent

Timeter® CAM4 Thermolectric Mist Tent Delivers Safety, Efficiency and Patient Comfort without CFCs Hoses connect the cooling unit to the tent, making it easier to adjust the height of the tent assembly, which is suspended over the patient bad.

- ▶ For use for children or adults
- ▶ Thermolectric cooling unit (no harmful CFC emissions)
- ▶ More efficient than refrigeration
- ▶ UL listed, CSA approved
- ▶ Adjustable tent height
- ▶ Sturdy base
- ▶ Quiet
- ▶ Incorporates the popular PCS 414 Compressor
- ▶ Feeder Tube Option
- ▶ Misto2gen® HV-12 nebulizer is capable of dispensing over 200cc of water per hour (3-1/3 cc/minute) to the patient.



10-T6393	Tent, Nebulizer, & Compressor
10-T5100	Tent, Nebulizer without Compressor
10-T4614	Compressor Only (Same as Above)
10-T2730	Nebulizer Only
10-T6308	Canopy Only