



# Vacuum/Pressure Gauges for CPAP or Neonatal Applications

▶ ± 30cm H<sub>2</sub>O in 1cm markings




**Bottom Mount**  
1/8" sm bore tubing  
Item #7-30-PBS




**Rear 1/4" NPT**  
Item #7-30-TRS

▶ ± 120cm H<sub>2</sub>O in 4cm markings




**Rear 1/4" NPT**  
Item #7-120-TRS


▶ ± 60cm H<sub>2</sub>O in 2cm markings



**Bottom Mount**  
1/8" sm bore tubing  
Item #7-60-PBS




**Rear 1/4" NPT**  
Item #7-60-TRS




**Bottom 1/4" NPT**  
Item #7-60-TBS

## NIF (Negative Inspiratory Force) Monitors with Adaptor

▶ ± 30cm H<sub>2</sub>O in 1cm markings




**Bottom Mount**  
1/8" sm bore tubing  
Item #7-30-PBR




**Rear 1/4" NPT**  
Item #7-30-TRR

▶ ± 120cm H<sub>2</sub>O in 4cm markings



**Rear 1/4" NPT**  
Item #7-120-TRR

▶ ± 60cm H<sub>2</sub>O in 2cm markings



**Bottom Mount**  
1/8" sm bore tubing  
Item #7-60-PBR



**Rear 1/4" NPT**  
Item #7-60-TRR




**Bottom 1/4" NPT**  
Item #7-60-TBR




**Gauge Cover**  
Clam Shell Design, Latex Free  
Item #7-NS-GC


All NIF Monitors are supplied with adaptors. Reusable replacement adaptors.



15mmID/22mmOD X 1/8" Small Bore Tubing Adaptor  
Item #7-30-A



1/4 NPT X 1/8" Small Bore Tubing Adaptor  
Item #7-00-A



15mmID/22mmOD X 1/4" NPT Tubing Adaptor  
Item #7-60-A