

# Nerve Stimulator and Accessories



## Compact Peripheral Nerve Stimulator


This compact peripheral nerve stimulator offers many functions found only in much larger units.

- ▶ Adjustable output voltage
- ▶ 1 pulse per 8 seconds
- ▶ 2 pulses per second
- ▶ Train of four (2 cycles of 2 per second)
- ▶ Tetanus frequency can be selected at 50 or 100 Hz
- ▶ Built-in battery test
- ▶ Accepts FDA compliant protected lead sets
- ▶ Uses standard 9 volt battery
- ▶ Designed for one-hand operation
- ▶ Visual indicator for pulse and battery condition

Size: 5.5" H x 1.5" W x .9" D    Warranty: 1 year  
Frequency: 1 Pulse/8 sec., 2 pulses/sec., 50 pulses/sec. (optional)  
Pulse Width: .30 millisecc    Output Current: 0 to 35 milliamp  
Battery Life (operating): up to 1000 hours  
Power: One standard 9v battery

### Replacement Parts:

Item #7-AA-57	54 in. Stimulator leads with alligator clips
Item #7-AA-67	54 in. Stimulator leads with snaps
Item #7-AA-10	54 in. Extension Set
Item #7-AA-62	Ball electrode set
Item #7-AA-09	Stimulator and Accessory Case




Item #7-AA-100    Peripheral Nerve Stimulator

Request a Sample

## Peripheral Nerve Stimulator Electrodes



Designed to be used for neuromuscular blockage monitoring. Provides for easy placement over the ulnar nerve on a patient's wrist. Dual element permits non-invasive monitoring using a single electrode site.



7-N6125101	1 pouch, 30 pouches/box, 1 box/case
7-N6125105	1/pouch, 30 pouches/box, 5 boxes/case
7-N612510	1/pouch, 30 pouches/box, 10 boxes/case